Antibody labeling kits

Contents

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Lumiprobe Corporation
(US and Worldwide)
9:00AM - 9:00PM EST

Phone: +1 888 973 63 53

Lumiprobe GmbH
(Germany and Europe)
8:00AM - 5:00PM CE(S)T

Phone: +49 511 16596811

Lumiprobe RUS Ltd
(Russia and CIS)
8:00 - 19:00 GMT+3 (Moscow)

Phone: +7 800 775 3271

ISO 9001:2015

Antibody labeling kits

The kit is used for the preparation of labeled antibodies having 2.5-3 labels per antibody. The kit contains all reagents and materials for ten labeling reactions, each for 100 µg of the antibody. The dye is a NHS ester, taken in controlled excess. It reacts with the amine groups of lysines of the antibody at mild alkaline pH. The purification of the antibody is achieved by gel filtration using spin columns. These kits contain sulfonated, water soluble dyes, that are most suitable for the labeling of sensitive proteins, including antibodies.

Kit components

| Kit component | Count | | | | |
|--|--------------|--------------|--------------|--------------|--------------|
| | 1321-10rxn | 3321-10rxn | 7321-10rxn | 5321-10rxn | 6321-10rxn |
| | 10 reactions | 10 reactions | 10 reactions | 10 reactions | 10 reactions |
| N1320, Sulfo-Cyanine3 NHS ester, 1 rxn | 10 | — | — | — | — |
| N3320, Sulfo-Cyanine5 NHS ester, 1 rxn | — | 10 | — | — | — |
| N7320, Sulfo-Cyanine5.5 NHS ester, 1 rxn | — | — | 10 | — | — |
| N5320, Sulfo-Cyanine7 NHS ester, 1 rxn | — | — | — | 10 | — |
| N6320, Sulfo-Cyanine7.5 NHS ester, 1 rxn | — | — | — | — | 10 |
| Desalting spin column, PBS | 10 | 10 | 10 | 10 | 10 |
| Desalting receptacle vial, 1.5 mL | 10 | 10 | 10 | 10 | 10 |
| Desalting waste vial, 2 mL | 10 | 10 | 10 | 10 | 10 |
| PBS tablet, for 100 mL of buffer | 1 | 1 | 1 | 1 | 1 |
| 15050, DMSO, labeling grade, 1 mL | 1 | 1 | 1 | 1 | 1 |
| 1584-05mL, Sodium azide solution, 3%, 0.5 mL | 1 | 1 | 1 | 1 | 1 |
| 1689-15mL, Sodium bicarbonate, 126 mg | 1 | 1 | 1 | 1 | 1 |

Store at temperature between +3 °C and +20 °C. Do not freeze!

Shelf life 12 months.

Protocol

1. Antibody preparation.

Optimal conditions for the labeling reaction are the following: antibody at a concentration of 1 mg/mL dissolved in 0.1 M sodium hydrocarbonate solution. Sodium azide is compatible with the labeling. If concentration of antibody is below 1 mg/mL, it should be concentrated. Antibody preparation should be free of aminoacids, proteins like BSA and buffer components that bring pH to unfavourable (2-7.5 or 9-12) ranges.

2. Setting up the reaction.

2.1. Add a solution of antibody in sodium hydrocarbonate (100 ug in 100 μ L*) to the fluorophore, vortex, and incubate for 30 min at room temperature.

* if a smaller amount of antibody is to be labeled, dissolve lyophilized dye in 10 μ L of anhydrous DMSO, and take 1 μ L of the solution per 10 ug of antibody. The solution of dye in DMSO should be consumed immediately.

3. Purification of the labeled antibody.


3.1. Prepare the column. The column should be at room temperature prior to use. The gel should be resuspended using vortex. Put the column into a waste receptacle (without cap), and centrifuge for 2 min at 1,000 g (speed must be carefully controlled — for standard 6 cm rotor, 1,000 g = 3800 rpm). The small bulge of the column should be oriented outwards.

3.2. Apply 400 μ L of 1x PBS buffer, centrifuge for 2 min at 1,000 g.

3.3. Remove the column from waste receptacle. Put it into collection tube (with cap). Apply 100 μ L of the reaction mixture onto the column, incubate for 1 min, and elute the antibody by centrifugation for 2 min at 1,000 g. Add 1/100 volume of 100x sodium azide. Antibodies can be aliquoted to improve their shelf life. Current aliquot can be stored at +4 °C the rest at -20 °C.

4. Measurement of dye to antibody ratio.

To calculate dye to antibody ratio, measure absorption spectrum of the conjugate, or absorption at 280 nm (A_{AB}), and at dye absorption maximum (A_{Dye}). A typical absorption spectrum of a dye labeled antibody is shown below. Depending on the dye, wavelength of the maximum may vary.

An absorption spectrum of a labeled antibody contains dye peak (with longer wavelength), and antibody peak (at around 280 nm). Dye to antibody ratio is calculated using the following formula:

$$\frac{\text{Dye}}{\text{AB}} = \frac{A_{Dye} \times \epsilon_{AB}}{(A_{AB} - A_{Dye} \times CF) \times \epsilon_{Dye}}$$

with **Dye/AB** — average number of fluorophores per antibody molecule, A_{dye} — optical density of the sample at dye absorption maximum, A_{AB} — sample optical density

at 280 nm, ϵ_{AB} — molar extinction coefficient of antibody at 280 nm (210,000 for IgG), ϵ_{Dye} — molar extinction coefficient of dye at absorption maximum (see table below), CF_{280} — correction factor for dye (see table below).

| Dye | λ_{max} nm | ϵ | CF_{280} |
|------------------|--------------------|------------|------------|
| sulfo-Cyanine3 | 548 | 162,000 | 0.06 |
| sulfo-Cyanine5 | 646 | 271,000 | 0.04 |
| sulfo-Cyanine5.5 | 673 | 195,000 | 0.11 |
| sulfo-Cyanine7 | 750 | 240,600 | 0.04 |
| sulfo-Cyanine7.5 | 778 | 222,000 | 0.09 |

Example of calculation

After labeling of IgG antibody with sulfo-Cyanine5 and purification, a solution has been obtained with absorption spectrum above. Determine dye to protein ratio.

Using absorption spectrum, the following values can be found: $A_{Dye} = 0.184$ at dye absorption maximum (646 nm, see the table), $A_{AB} = 0.066$ (at 280 nm), $\epsilon_{AB} = 210000$ for IgG; $\epsilon_{Dye} = 271,000$, and $CF_{280} = 0.04$.

$$\frac{Dye}{AB} = \frac{A_{Dye} \times \epsilon_{AB}}{(A_{AB} - A_{Dye} \times CF) \times \epsilon_{Dye}} = \frac{0.184 \times 210\,000}{(0.066 - 0.184 \times 0.04) \times 271\,000} = 2.43$$

Dye/AB — is 2.43 fluorophore molecules per antibody.

Interpretation of results, and antibody storage

In most cases, optimal dye to antibody ratio is 2-3 dye molecules per antibody. Increase of dye loading above this value does not improve fluorescent signal because of the concentrational quenching of the fluorescence. When too few fluorophores are attached during the reaction, decrease antibody loading per reaction. Use of expired kit can also lead to this result.

We recommend to test binding of labeled antibodies with their antigens. The labeled antibodies can be stored at $-20\text{ }^{\circ}\text{C}$. An aliquot that is currently in use, should be stored at $+4\text{ }^{\circ}\text{C}$ to avoid freeze-thaw cycles. The stability of the conjugate is determined by antibody itself rather than fluorophore or the chemical bond. Labeled antibodies should not be exposed to direct sunlight for extended time, but tolerate ambient lighting well.

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Dye/AB = 2,43 Fluorophormoleküle pro Antikörpermolekül.
Issued by TÜV NORD CERT GmbH