PKH Cell Membrane Labeling Kits manual

Lumiprobe Corporation. All rights reserved.

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Shelf life 12 months.

Recommendations for using the kit

- Optimal concentrations of the dye and cells can vary depending on cell and study type, so evaluate cell viability, homogeneity, and fluorescence intensity after staining.
- Do not use azide-containing solutions when staining with PKH dyes.
- Staining is more homogeneous when cell suspension is used.

Protocol

Protocol for cell membrane labeling with PKH dyes for RAW264.7 adhesion culture, 1×10^6 cells/sample, the final concentration of PKH dye 2 μM , final volume 200 μL .

1. Prepare PKH dye solution immediately before staining. Add 1 μL of PKH dye solution (*PKH dye, 1 mM solution in isopropanol*) to 9 μL of 96% ethanol, and add 4 μL of the resulting solution to 100 μL of *PKH Dyes Diluent, 1x*.

**1x PKH Dyes Diluent is available either in ready-to-use form (K6150, PKH Dyes Diluent, 1x) or as a 5x concentrate (K7150, PKH Dyes Diluent, 5x). To dilute 5x PKH Dyes Diluent, use sterile bidistilled water.*
2. Remove the cell culture from the surface with a scraper in Hanks' solution (HBSS). Count the cells in the sample. Add 3 mL of Hanks' solution, centrifuge at 400 $\times g$ for 6 min at room temperature.

**Serum proteins and lipids also bind the dye, so it is recommended to wash the cells once with serum-free medium or phosphate buffer saline.*
3. Remove the supernatant with a pipette, and resuspend a necessary amount of cells (e. g. 1×10^6 cells) in 100 μL of *PKH Dyes Diluent, 1x*. Add 100 μL of PKH dye solution prepared in step 1. Pipette and allow to stand at room temperature for 5 min. The final PKH dye concentration in the cell solution is 2 μM .

**To obtain reproducible results, minimize the volume of the supernatant before cell resuspending.*

**Do not leave cells in PKH Dyes Diluent for a long period.*

**Staining is almost instant, so rapid cell dispersion in the dye solution is important to produce bright homogeneous and reproducible labeling.*
4. Add 2 mL of fetal bovine serum to stop the reaction, and incubate for 1 min. Centrifuge at 400 $\times g$ for 10 min at room temperature.

**To stop the reaction, do not use a serum-free medium or buffered salt solution that results in dye aggregates.*

5. Remove the supernatant, resuspend the cells in 5 mL of complete culture medium, and transfer them to a new tube. Take aliquots to evaluate cell viability with trypan blue. Centrifuge at 400 x g for 10 min at room temperature.
6. Resuspend the cells in a buffer for further analysis (microscopy, flow cytometry, etc.).

**Stained cells can be fixed with 2% paraformaldehyde, and staining remains stable for at least 3 weeks if samples are protected from light.*

Recommendations for storage

PKH dye solution can be stored at room temperature or in a refrigerator protected from light. Check the solution for precipitation before use. If a precipitate is seen in the dye solution, slightly warm it in a water bath at 37 °C and ultrasonicate or vortex until redissolved.

PKH Dyes Diluent is delivered as a 1x or a 5x solution in a sterile container. Store in a refrigerator and adjust to room temperature immediately before use.

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