## LumiPure<sup>™</sup> Plasmid spin miniprep Kit manual

## Contents

English: LumiPure Plasmid spin miniprep Kit manual .....	3-7
Deutsch: Handbuch für das LumiPure spin Miniprep-Kit für Plasmide	8-12
Русский: Инструкция к набору LumiPure для выделения плазмидной ДНК .....	13-17

# LumiPure Plasmid spin miniprep Kit manual

The kit allows rapid (about 15 minutes) and highly efficient purification of plasmid DNA (up to 30 µg) from *Escherichia coli* cultures using spin columns. The purified DNA is suitable for any subsequent application, such as PCR, restriction enzyme digestion, ligation and transformation, sample preparation for Sanger sequencing and NGS, etc.

## Kit components

Kit component	Count
	<b>21583</b>
	<b>50 minipreps</b>
M1250, Resuspension Solution, 15 mL	1
31650, RNase A (10 mg/mL), 150 µL	1
M1550, Lysis Solution, 15 mL	1
M1450, Neutralization Buffer, 20 mL	1
P2450, Wash Solution A (with GuHCl), 30 mL	1
K2250, Wash Solution B (Concentrate to be diluted 5x with ethanol 96%), 10.0 mL	2
K1350, Elution Buffer (10 mM Tris-HCl, pH 8.5), 10 mL	1
Spin column CNP30-A (up to 30 µg, blue ring)	50
Collection tube for spin column, 2 mL	50

Store at room temperature.

Shelf life 12 months.

## Equipment and reagents supplied by user:

- microcentrifuge with rotor for 1.5 mL tubes and capable of centrifuging  $>10,000$  RPM ( $6,700 \times g$ );
- 96 % ethanol;
- 1.5 mL microcentrifuge tubes (2 tubes for DNA purification from 1 biological sample);
- (optional) tabletop centrifuge with rotor for 15 mL tubes and capable of centrifuging  $>4,500 \times g$  for centrifugation of overnight bacterial culture (alternatively, a microcentrifuge with rotor for 1.5 mL tubes can be used).

## Before starting

1. Transfer the contents of *RNase A* vial into the *Resuspension Solution* vial and mix thoroughly. Mark on the label that *RNase A* has been added.
2. Dilute the concentrate of *Wash Solution B* 5x with 96 % ethanol (add 4 volumes of 96 % ethanol to 1 volume of concentrate specified on the bottle). Mark on the label that ethanol has been added.
3. If any precipitate has formed in *Lysis Solution*, *Neutralization Buffer*, *Wash Solution A*, incubate the solutions at a temperature not higher than 50 °C until the precipitates have been completely dissolved. Before use let the solutions cool down to 25 °C.

## Purification of plasmid DNA

All centrifugation steps are carried out at room temperature, at 10,000–13,000 RPM (6,700–11,000 x g) (unless specified otherwise).

Process bacterial cultures with OD<sub>600</sub> no more than 30 AU (for instance, 15 mL of the culture with OD<sub>600</sub> = 2 AU). If a larger cell amount should be processed, volumes of *Resuspension Solution*, *Lysis Solution*, *Neutralization Buffer* should be increased proportionally, and the supernatant should be loaded to the column in several stages.

1. Centrifuge 2.5–7 mL of *E. coli* overnight culture (10–15 mL for low-copy-number plasmids) for 5 minutes, 5,000 RPM (4,500 x g) or for 1 minute, 13,000 RPM (10,000 x g). Remove the supernatant and all remaining culture medium.
2. Resuspend the cell pellet in 250 µL of *Resuspension Solution*. Transfer the cell suspension into a clean 1.5 mL microcentrifuge tube.
3. Add 250 µL of *Lysis Solution*, mix gently by inverting the tube 4–6 times. Check that a turbid solution has become viscous and clear, and immediately proceed to the next step.

*! In case of processing 10–15 mL of bacterial culture, mix thoroughly by inverting the tube until the solution has become viscous and clear. In some instances, more inversions may be required.*

*! Do not vortex the sample. This can cause shearing of chromosomal DNA.*

4. Add 350 µL of *Neutralization Buffer* to the sample, mix gently by inverting the tube 4–6 times. As a result precipitations in the form of flakes should appear in the sample.

*! In case of processing 10–15 mL of bacterial culture, after the tube contents has been mixed, shake the tube for a short period of time (a few seconds).*

5. Centrifuge the sample for 5 minutes.
6. Place a spin column in a collection tube. Carefully transfer the supernatant to the column. Centrifuge the column for 30 seconds. Discard the flow-through and place the column back into the same collection tube.

*! (optional) To remove endotoxins or trace nuclease activity (when working with EndA+ strains) add 500  $\mu$ L of Wash Solution A to the column (this step can result in decreasing of DNA yield by up to 20 %). Centrifuge the column for 30 seconds. Discard the flow-through and place the column back into the same collection tube.*

7. Add 500  $\mu$ L of *Wash Solution B*, centrifuge the column for 30 seconds. Discard the flow-through and place the column back into the same collection tube.
8. Add 500  $\mu$ L of *Wash Solution B*, centrifuge the column for 3 minutes. Discard the collection tube with the flow-through.
9. Place the spin column into a clean 1.5 mL microcentrifuge tube. Add 50–100  $\mu$ L of *Elution Buffer* to the center of the spin column membrane. Incubate at room temperature for 1 minute. Centrifuge the column for 2 minutes. The microcentrifuge tube contains the purified DNA.

*! For increased DNA concentrations, use a low volume of Elution Buffer. Elution with volumes of less than 50  $\mu$ L is not recommended because this can be insufficient to entirely wet the membrane resulting in low DNA recovery. For increased DNA yield, use a higher volume of Elution Buffer (100  $\mu$ L).*

*! If necessary, DNA can be eluted with deionized water.*

## Note

When measuring DNA concentration, dilute the sample only with TE buffer pH 8.5 or *Elution Buffer* supplied with the kit. Otherwise, DNA concentration and purity ( $OD_{260}$  and  $OD_{260}/OD_{280}$ , respectively) can be estimated incorrectly.

**Storage of purified DNA:** for long-term storage  $-20\text{ }^{\circ}\text{C}$ , for short-term storage  $+4\text{ }^{\circ}\text{C}$ .

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## **Lumiprobe Corporation**
## **Lumiprobe GmbH**
## **Lumiprobe Limited**
Issued by TÜV NORD CERT GmbH