

Staining Live Cells with Endoplasmic Reticulum Stains LumiTracker ER

LumiTracker ER Green and LumiTracker ER Red are cell-permeant, highly selective for the endoplasmic reticulum (ER) stains that can be used for live cell imaging. Both ER stains are derivatives of BDP dyes coupled with Glibenclamide (Glyburide). Glibenclamide binds to the sulphonylurea (SUR) receptors of ATP-sensitive potassium channels (K_{ATP}), which are prominent on the endoplasmic reticulum.

The staining is partially retained after fixation with formaldehyde. LumiTracker ER Green and LumiTracker ER Red are unsuitable for staining cells after fixation.

Note: The pharmacological activity of glibenclamide could potentially affect ER function. Variable expression of sulphonylurea receptors in some specialized cell types may result in non-ER labeling.

Preparation of solutions

1.1 Preparation of the stock solutions

- Dissolve 50 μ g LumiTracker ER Green in 64 μ L DMSO to obtain 1 mM of stock solution.
- Dissolve 50 μ g LumiTracker ER Red in 55 μ L DMSO to obtain 1 mM of stock solution.

Note: We recommend storing the stock solutions at $-20\text{ }^{\circ}\text{C}$ or $-80\text{ }^{\circ}\text{C}$ away from light. Avoid repeated freeze-thaw cycles.

1.2 Preparation of working solution

Dilute the stock solution of LumiTracker ER stain in Hank's Balanced Salt Solution with calcium and magnesium (HBSS/Ca/Mg) to obtain the working solution. Use the concentration from 100 nM to 1 μ M. The dilution of the LumiTracker ER stain depends on the cell type and density and should be defined experimentally.

Cell staining

2.1 Suspension cell staining

1. Centrifuge cell suspension at 1000 g for 3-5 min; discard the supernatant.
2. Wash cells twice with prewarmed HBSS, 5 min each time at $37\text{ }^{\circ}\text{C}$.
3. The recommended cell density is 1×10^6 cells/mL.
4. Add 1 mL of prewarmed LumiTracker ER working solution to the tube, and incubate cells for 15-30 min at $37\text{ }^{\circ}\text{C}$ and 5% CO_2 .
5. Centrifuge cells at 400 g for 3-4 min at room temperature; discard the supernatant.
6. Wash stained cells twice with medium.
7. Resuspend cells with serum-free cell culture medium or PBS.

- Analyze cells by flow cytometry.
- For fluorescent microscopy, transfer a drop of the stained cell suspension onto a glass slide. Cover the cells with a coverslip.

2.2 Adherent cell staining

- Grow cells on a sterile coverslip.
- Adherent cells can be stained directly on the coverslip.
- Aspirate the medium from coverslip.
- Rinse cells with prewarmed HBSS at 37 °C.
- Add 100 µL of prewarmed LumiTracker ER working solution and gently shake it to cover the cells completely. Incubate cells for 15-30 min at 37 °C and 5% CO₂.
- Wash stained cells twice with medium.
- If detected by flow cytometry, cells need to be resuspended before staining.
- For fluorescent microscopy, invert the coverslip with stained cells onto the slide to place them between the slide and the coverslip.

Cell fixation

- If stained cells are to be fixed, use an incubation in 4% formaldehyde for 2 min at 4 °C.
- Wash fixed cells twice in PBS for 5 min each time.
- Mount cells under a coverslip using a [mounting medium](#).

Lumiprobe Corporation

201 International Circle, Suite 135
Hunt Valley, Maryland 21030
USA
Phone: +1 888 973 6353
Fax: +1 888 973 6354
Email: order@lumiprobe.com

Lumiprobe GmbH

Feodor-Lynen-Strasse 23
30625 Hannover
Germany
Phone: +49 511 16596811
Fax: +49 511 16596815
Email: de@lumiprobe.com

Lumiprobe RUS Ltd

Kotsyubinsky street, 4
121351 Moscow
Russian Federation
Phone: +7 800 775 3271
Email: ru@lumiprobe.com