LumiCell PKH26 cell membrane labeling kit manual

Contents

Deutsch: Handbuch für das LumiCell PKH26 Kit für die Zellmembranmarkierung	3-5
Русский: Инструкция к набору для мечения клеточных мембран LumiCell PKH26	6-8
English: LumiCell PKH26 cell membrane labeling kit manual	9-11

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Shelf life 12 months.

Recommendations for using the kit

- Optimal concentrations of the dye and cells can vary depending on cell and study type, so evaluate cell viability, homogeneity and fluorescence intensity after staining.
- Do not use azide-containing solutions when staining with PKH26 dye.
- Staining is more homogeneous when cell suspension is used.

Protocol

Protocol for cell membrane labeling with PKH26 for RAW264.7 adhesion culture, 1×10^6 cells/sample, the final concentration of PKH26 $2 \mu\text{M}$, final volume $200 \mu\text{L}$.

1. Prepare PKH26 solution immediately before staining. Add $1 \mu\text{L}$ of 1mM PKH26 solution to $9 \mu\text{L}$ of 96% ethanol, add $4 \mu\text{L}$ of the resulting solution to $100 \mu\text{L}$ of PKH26 Diluent, 1x.
2. Remove the cell culture from the surface with a scraper in Hanks' solution (HBSS). Count the cells in the sample. Add 3mL of Hanks' solution, centrifuge at 400g for 6min at room temperature.

**Serum proteins and lipids also bind the dye, so it is recommended to wash the cells once with serum-free medium or phosphate buffer saline.*

3. Remove the supernatant with a pipette, resuspend a necessary amount of cells (e. g. 10^6 cells) in $100 \mu\text{L}$ of PKH26 Diluent, 1x. Add $100 \mu\text{L}$ of PKH26 solution prepared during step 1. Pipette and allow to stand at room temperature for 5min . The final PKH26 concentration in the cell solution is $2 \mu\text{M}$.

**To obtain reproducible results, minimize the volume of the supernatant before cell resuspending.*

**Do not leave cells in PKH26 Diluent for a long period.*

**Staining is almost instant, so rapid cell dispersion in the dye solution is important to produce bright homogeneous and reproducible labeling.*

4. Add 2mL of fetal bovine serum to stop the reaction, incubate for 1min . Centrifuge at 400g for 40min at room temperature.

**To stop the reaction, do not use a serum-free medium or buffered salt solution that results in dye aggregates.*

5. Remove the supernatant, resuspend the cells in 5mL of complete culture medium, transfer to a new tube. Take aliquots to evaluate cell viability with trypan blue. Centrifuge at 400g for 40min at room temperature.
6. Resuspend the cells in a buffer for further analysis (microscopy, flow cytometry


etc.).


**Stained cells can be fixed with 2% paraformaldehyde, and staining remains stable for at least 3 weeks if samples are protected from light.*

Recommendations for storage

PKH26 solution can be stored at room temperature or in a refrigerator protected from light. Check the solution for precipitation before use. If a precipitate is seen in the dye solution, slightly warm it up on a water bath at 37°C and ultrasonicate or vortex until redissolved.

PKH26 Diluent is delivered as a 1x or a 5x solution in a sterile container. Store in a refrigerator and adjust to room temperature immediately before use.

Lumiprobe Corporation

(US and Worldwide)

9:00AM - 9:00PM EST

201 International Circle, Suite 135

Hunt Valley, Maryland 21030

USA

Phone: +1 888 973 63 53

Lumiprobe GmbH

(Germany and Europe)

8:00AM - 5:00PM CE(S)T

Feodor-Lynen-Strasse 23

30625 Hannover

Germany

Phone: +49 511 16596811

Lumiprobe RUS Ltd

(Russia and CIS)

8:00 - 19:00 GMT+3 (Moscow)

Kotsyubinskogo 4, bld. 3

121351 Moscow

Russian Federation

Phone: +7 800 775 3271

ISO 9001:2015

The certificate No.: 44 100 202070

Issued by TÜV NORD CERT GmbH

www.lumiprobe.com