

Magnetic Beads (100 mg/mL)

<http://www.lumiprobe.com/p/magnetic-beads>

Magnetic silica-based particles represent the modern standard for rapid and efficient DNA and RNA isolation. This technology delivers preparations of the highest purity, free from inhibitors and with minimal sample loss. Offering seamless compatibility with automated platforms, it is equally suited for manual protocols and high-throughput research. You will obtain pure, ready-to-use nucleic acids, ideal for even the most demanding downstream applications.

General properties

Legal statement: This Product is offered and sold for research purposes only. It has not been tested for safety and efficacy in food, drug, medical device, cosmetic, commercial or any other use. Supply does not express or imply authorization to use for any other purpose, including, without limitation, in vitro diagnostic purposes, in the manufacture of food or pharmaceutical products, in medical devices or in cosmetic products.